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Preventive maintenance schedule

Instrument	Model & Manuf.	Action #	Action	Frequency
Flame Photometer	Corning 410C	FP1	Empty water container	Upon each use
			Check air line for condensation	Upon each use
			Check U tubes filled with deionized water	Upon each use
		FP2	Check nebulizer	Weekly
		FP3	Check: 1. Head drain 2. Mixing chamber 3. U tube and drain tube and clean if necessary 4. Refill the U tube with Deionized water	Monthly
		FP4	1. Check air and fuel tubing and connection for leaks.	Every 6 months
Ion Chromatography	Dionex ICS-2100	ICD1 ICD2	1. Calibration 2. Replacement of eluent generator	Every 6 months Annually
Atomic Absorption Spectroscopy/Graphite Furnace	Shimadzu AA680	GFSHIM	1. Clean graphite cap, graphite holder, window plates, and emission detector hole and temperature sensor unit.	Monthly
Ratio Turbidimeter	HACH	TURH	1. Calibration.	Every 6 months
Turbidimeter	WTW	TURW	1. Calibration.	Every 6 months
Conductivity Meter	Jeneway 4320	COND	1. Calibration.	Annually
pH-Meter	WWW Inolab	PH1	1. Replace the battery	When needed
		PH2	1. Calibration.	Weekly
UV Spectrophotometer	Shimadzu 160 A	UVSH1	1. Calibration.	Every 3months
Ultra pure water unit	MILLIPORE	UM1 UM2	1. Check the conductivity and PH 2. Check the filters and UV.	Every day When needed
Electronic balance1	Shimadzu AX200	EBAL	1. Dustening, and clean the balance using neutral detergent 2. Avoid using organics solvents, chemicals, or dusting sprays	Upon each use
All laboratory instruments			Clean dust and visual inspection	weekly
Water treatment system	CHiyi	WT 1 WT 2 WT3	Change the filter kit. Check Salt amount and addition of salt. Check the system by measuring the PH and conductivity.	Every 6 months. Every week. Every day (Data recorded on Bench Sheet for TDS,PH Daily)
UV.Spectrophotometer	Merk prove 300	UV 300-1	Calibration.	Every 6 months.

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Air density plates: Every three weeks

Swab method : Every week


Distilled Water Sterility: Every Month

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Refrigerator No.:

Correction Factor:

Date	Temp. after correction C°	Time	Name & Signature

LAB5

MASTER COPY